


HYGIENE & DISINFECTION PRODUCTS


— Reliable Disinfection —



☒ FUNGICIDAL ☒ BACTERICIDAL ☒ TUBERCULOCIDAL ☒ VIRUCIDAL ☒ MYCOBACTERICIDAL ☒ YEASTICIDAL



The dynamic, forward looking, profession team of Prevest DenPro Limited Research & Development department has successfully developed a comprehensive range of revolutionary latest generation of broad spectrum, rapid action cleaning & disinfection solutions for dental and medical equipments. The products are recommended for use for cleaning and disinfection of all types of non invasive medical devices used in dental and medical operatories. The products have been tested for the bactericidal, virucidal, fungicidal, Sporicidal and yeasticidal efficacy. The products have been formulated with latest generation of patented blend of quaternary ammonium compounds, amines, surfactants and excipients which are highly effective against gram negative and gram positive bacteria, antibiotic resistant bacteria, human virus, veterinary viruses, fungi, tuberculosis. The products are free of phenol, aldehydes, halogens and chlorine. The product formulations are highly effective, easy to use and safe.



ViroGuard

Surface Disinfectant Spray for Medical Devices

VIROGUARD is a rapid-action, broad-spectrum antimicrobial spray formulated specifically for cleaning and disinfecting the surfaces of medical equipment. Regular application between patient consultations helps minimize the risk of cross-contamination in medical and dental settings. It is ideal for disinfecting high-contact surfaces such as benches, examination tables, dental units, handles, spittoons, and more. VIROGUARD offers quick action with powerful cleaning and disinfecting performance, complemented by a mild, refreshing scent for comfortable use.

Key features

- Bactericidal, fungicidal, tuberculocidal, mycobactericidal, virucidal
- Does not contain phenols, aldehydes, chlorine or EDTA
- Does not affect the medical equipment
- Fast acting
- Leaves no residue after drying

ACTIVITY SPECTRUM	STANDARD	STRAINS
BACTERICIDAL (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae
FUNGICIDAL (Dirty conditions)	EN 13624	Candida albicans



Recommended usage:

Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 sec for effective disinfection and wipe with the clean lint free disposable cloth. (Note : Do not rinse)

Indications for use:

It is used for cleaning, decontaminating and disinfecting surfaces (such as trolleys, beds, chairs, trays, work benches, examination tables, dental units, handles, spittoons, other non invasive equipments and instruments).

Compatibility

Due to its low alcohol content, Viroguard is friendly to sensitive surfaces while simultaneously providing extremely rapid disinfection times.

Composition

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant , IPA, pH modifier and water

Packaging

[REF] 16012: 500 ml Spray Bottle,

[REF] 16013: 1 Litre Can

[REF] 16014: 5 Litre Can

Physical Properties	
Appearance	Transparent Solution
Density	0.9 -1.01g/cm³ at 20°C
pH	9.0 - 10 at 20°C
Odour	Mild (alcohol)
Storage	5°C - 35°C
Stability	3 Years



ViroShield

Disinfectant for Suction System

VIROSHIELD is a powerful concentrated formulation intended for the thorough cleaning and disinfection of common vacuum lines, amalgam filters, and spittoons. It effectively clears organic matter such as saliva, blood, protein deposits, and debris from drilling procedures, helping to reduce biofilm buildup and ensuring hygienic internal pipe surfaces. The solution leaves behind a subtle, pleasant fragrance while eliminating foul odors. Its wide-ranging antimicrobial properties and strong cleansing ability supports a hygienic and safe setting for both healthcare professionals and patients. For best hygiene practices, it is advised to use VIROSHIELD daily at the end of each workday.

Key features

- Cost effective-2% dilution
- Bactericidal, fungicidal and Virucidal
- Eradicates bad smell
- Does not foam
- Does not contain phenols, aldehydes or chlorine
- Contact Time: 15 Minutes

ACTIVITY SPECTRUM	STANDARD	STRAINS	CONTACT TIME
BACTERICIDAL (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	5 Min.
FUNGICIDAL (Dirty conditions)	EN 13624	Candida albicans Aspergillus brasiliensis	5 Min. 60 Min.



Recommended usage:

The concentrate is to be diluted with water to 2% concentration. Pass the solution through the suction system once or twice a day. Leave for 15 minutes to take the effect. After disinfection, rinse the suction system with the portable water.

Indications for use:

It is indicated for daily hygiene and disinfection of dental suction systems as a part of their general disinfecting routine.

Compatibility

Viroshield is compatible with all common types of vacuum systems, amalgam separators and spittoons.

Composition

Didecyl-dimethyl ammonium chloride, N-(3-aminopropyl)-N-dodecylpropano-1,3-Diamine, Quaternary Ammonium Salts, Tetra Potassium Pyrophosphate, Isotridecanol Ethoxylated, 2-Aminoethanol, EDTA, Potassium Hydroxide, Dipropylene Glycol Methyl Ether and IPA.

Packaging

[REF] 16001: 1 Litre Dosing Bottle,

[REF] 16002: 5 Litre Can

Physical Properties

Appearance	Pale Yellow Solution
Density	0.9 - 1.01g/cm³ at 20°C
pH (2%w/v)	9.0 - 10.0 (Neutral at 20°C)
Odour	Mild Aromatic (Lemon)
Storage	5°C - 35°C
Stability	3 Years



ViroZym

Enzymatic Concentrate for Cleaning and Disinfection of Medical Devices

VIROZYM is a high-strength enzymatic cleaner specifically formulated for the pre-disinfection or pre-sterilization cleaning of surgical instruments and endoscopes. Its potent composition include a protease that cleans peptide bond. This enhances the removal of most resistant organic contaminants. Gentle on the surfaces, VIROZYM offers a safe cleaning process for users. It is suitable for use in both ultrasonic cleaners and immersion baths.

Key features

- Effective with 1% dilution.
- Bactericidal, fungicidal.
- Cleans, disinfects and deodorizes in one step.
- Removes stubborn deposits.
- Does not foam
- Does not contain phenols, aldehydes or chlorine

ACTIVITY SPECTRUM	STANDARD	STRAINS	CONTACT TIME
BACTERICIDAL* (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	5 Min.
FUNGICIDAL (Dirty conditions)	EN 13624	Candida albicans	5 Min.



Recommended usage:

The concentrate is to be diluted with portable water to 1% concentration. Immerse instruments immediately after use in the solution and soak for 15 minutes. Rinse the instruments with water and dry after disinfection. Discard the used solution daily and prepare fresh solution daily. Do not mix with other cleaning agents or disinfectants.

Indications for use:

It is used for manual cleaning and disinfection of dental and medical devices such as rotary and surgical instruments respectively. It can also be used in ultrasonic cleaners for high level enzymatic cleaning and disinfection of dental and medical instruments.

Compatibility

Virozym is compatible with most materials such as stainless steel, aluminium, glass, ceramics, hard plastics, ebonite etc.

Composition

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Poly(oxy-1,2-ethanediyl), alpha-[2-(didecylmethylammonio) ethyl]-omegahydroxy-propanoate, Didecylidimethylammonium chloride, Isotridecanol, ethoxylated, 2-(2-Butoxyethanol), Subtilisin, thickening agents, emulsifiers and fragrance.

Packaging

REF 16005: 1 Litre Dosing Bottle

REF 16006: 5 Litre Can

Physical Properties

Appearance	Light Yellow Solution
Density	0.9 - 1.01g/cm ³ at 20°C
pH (1% w/w)	9.0 - 10.0 (Neutral at 20°C)
Odour	Neutral
Storage	5°C - 35°C
Stability	3 Years



ViroGuard Foam

Surface Disinfectant Foam Spray for Medical Devices

Viroguard Foam is a foaming disinfectant spray with a broad antimicrobial spectrum of activity. It combines excellent cleaning and disinfecting properties and ensures the protection of patients and staff. It's alcohol-free formulation is ideal for the rapid disinfection of sensitive medical equipment surfaces, as incubators, monitors, plexiglass, upholstery materials. It is also suitable for all surfaces of medical equipment such as beds, stretchers, work benches and all other medical devices.

Key features

- Fast acting, ready-to-use alcohol free spray foam disinfectant
- Bactericidal, fungicidal
- Leaves no residue after drying
- Does not affect medical equipments
- Suitable for alcohol sensitive materials such as leatherite, plastic
- Does not contain phenol, aldehydes and chlorine
- Excellent wetting properties

ACTIVITY SPECTRUM	STANDARD	STRAINS
BACTERICIDAL* (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae
FUNGICIDAL (Dirty conditions)	EN 13624	Candida albicans



Recommended usage:

Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 sec for effective disinfection and wipe with the clean lint free disposable cloth. For difficult surfaces such as keyboards, displays, monitors and control panels, use wet clean disposable cloth with spray foam and thoroughly wipe the surfaces.

Indications for use:

- It is an alcohol free foam spray used for disinfection of alcohol sensitive non invasive medical devices such as dental units, beds, stretchers, work benches etc.
- It can also be used for disinfection of non invasive medical equipments such as monitors, displays, keyboards and control panels etc.

Compatibility

Viroguard Foam is compatible with most materials such as stainless steel, aluminium, glass, ceramics, hard plastics, ebonite, plexiglass, upholstery materials.

Composition

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant, Sodium lauryl sarcosinate, pH modifier and water


Packaging

[REF] 16009: 500 ml Spray Bottle,

[REF] 16010: 1 Litre Can

[REF] 16011: 5 Litre Can

Physical Properties	
Appearance	Foaming Solution
Density	0.9 - 1.01g/cm ³ at 20°C
pH	9.0 - 10.0 at 20°C
Odour	Neutral
Storage	5°C - 35°C
Stability	3 Years



ViroWipes

Rapid Disinfection Wipes for Medical Devices

ViroWipes is a one-step, ready to use, disposable wipes for the rapid disinfection of all kinds of wipeable surfaces of non-invasive medical devices such as hand pieces, contra-angle, turbines, scaler inserts and other non invasive instruments and equipments with increased risk of infection.

Key features

- One-step-simply wipe to clean and disinfect surfaces of non-invasive medical devices.
- Revolutionary formula-impregnated with patented blend of antimicrobial compounds in alcoholic solution.
- Leaves no residue on the surfaces
- Pleasant fragrance

ACTIVITY SPECTRUM	STANDARD	STRAINS
BACTERICIDAL* (Dirty conditions)	EN 16615	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae
FUNGICIDAL (Dirty conditions)	EN 16615	Candida albicans



Recommended usage:

Wipe the equipment with ViroWipes and massage thoroughly over the surface and provide scrubbing action to ensure good surface contact with the disinfectant. Leave the instruments/equipments to dry. Do not rinse with water and do not mix with any other disinfectant.

Indications for use:

It is used for a rapid disinfection and cleaning of medical devices and surfaces like hand pieces, contra-angles, turbines, scaler inserts etc.

Compatibility

ViroWipes is compatible with most materials such as dental chairs, trays, handpiece surfaces, light handles, countertops, X-ray equipment, and other non-invasive medical device surfaces.

Composition

N-(3-aminopropyl)-N-dodecylpropano-1,3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant , IPA, pH modifier, water

Packaging

REF 16015: 80 x wipes mono pack

Physical Properties

Appearance	Colored Wipes
pH	9.0 - 10.0 at 20°C
Odour	Neutral
Storage	5°C - 35°C
Stability	3 Years

HYGIENE & DISINFECTION PRODUCTS

Reliable Disinfection



Prevest DenPro Limited

Export Promotion Industrial Park,
Bari Brahmana, Jammu - 181133 (India)

+91 88990 74151 (Customer Care Number)
+91 88993 36171 (For International Market)
+91 98585 11041 (For Indian Market)

info@prevestdenpro.com
customercare@prevestdenpro.com

www.prevestdenpro.com